


ARBEITSSCHEMA



SCHNEIDAUFSATZ Ohne Rotationsschutz

Plattform	Höhe (mm)	Artikelnr.
Narrow	3.5	IPD/KA-IN-35
	5.0	IPD/KA-IN-50
	6.5	IPD/KA-IN-65
Standard	3.5	IPD/KA-IR-35
	5.0	IPD/KA-IR-50
	6.5	IPD/KA-IR-65
Wide	3.5	IPD/KA-IW-35
	5.0	IPD/KA-IW-50
	6.5	IPD/KA-IW-65

Version mit Fenster für gewinkelten Schraubenkanal nach Verfügbarkeit.



Mit Rotationsschutz

Plattform	Höhe (mm)	Artikelnr.
Narrow	3.5	IPD/KR-IN-35
	5.0	IPD/KR-IN-50
	6.5	IPD/KR-IN-65
Standard	3.5	IPD/KR-IR-35
	5.0	IPD/KR-IR-50
	6.5	IPD/KR-IR-65
Wide	3.5	IPD/KR-IW-35
	5.0	IPD/KR-IW-50
	6.5	IPD/KR-IW-65

Version mit Fenster für gewinkelten Schraubenkanal nach Verfügbarkeit.



BEFESTIGUNGSSCHRAUBE


	Schlüssel	Artikelnr.
Alle Schneidaufsätze	Sechskant 1.20	IPD/KA-TI-00



SCHNEIDAUFSATZ-SETS

Plattform	Heights (mm)		
Narrow	3.5 / 5.0 / 6.5	IPD/KT-IN-00	IPD/KT-IN-01
Standard	3.5 / 5.0 / 6.5	IPD/KT-IR-00	IPD/KT-IR-01
Wide	3.5 / 5.0 / 6.5	IPD/KT-IW-00	IPD/KT-IW-01

Alle Sets inklusive Befestigungsschraube IPD/KA-TI-00.



SCHLÜSSEL SCAN BODY

	Höhe Body height (mm)	Artikelnr.
	10	IPD/PH-20-18
	15	IPD/PH-20-25



ipd ist eine eingetragene Marke von IMPLANT PROTESIS DENTAL 2004 SL. Die übrigen genannten Marken sind Eigentum ihrer jeweiligen Eigentümer, wobei IMPLANT PROTESIS DENTAL 2004 SL nicht mit ihnen oder den genannten Marken in Verbindung steht.

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Nobel Biocare® Branemark System® Klockner® KL™	Ø3.5	M. Rot.Schutz	1.5	IPD/AA-IN-02/3D	IPD/KT-IN-00
		O. Rot.Schutz		IPD/AA-IN-03/3D	IPD/KT-IN-01
Nobel Biocare® Branemark System® Biomet® 3i® Osseotite® BTI® External Hex Sweden y Martina® Outlink® Microdent® Microdent® Universal™ Klockner® KL™ Neodent® Helix® HE	Ø4.1	M. Rot.Schutz	0.5	IPD/AA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/AA-IR-02/3D	
			2.5	IPD/AA-IR-04/3D	
			0.5	IPD/AA-IR-01/3D	
	S3-S4	O. Rot.Schutz	1.5	IPD/AA-IR-03/3D	IPD/KT-IR-01
			2.5	IPD/AA-IR-05/3D	
Nobel Biocare® Branemark System®	Ø5.1	M. Rot.Schutz	1.3	IPD/AA-IW-02/3D	IPD/KT-IW-00
		O. Rot.Schutz		IPD/AA-IW-03/3D	-
Nobel Biocare® Multi-unit System	Ø4.8	M. Rot.Schutz	0.5	IPD/AB-IR-00/3D	-
			0.5	IPD/AB-IR-02/3D	
		O. Rot.Schutz	0.2	IPD/AB-IR-01/3D	-
			0.2	IPD/AB-IR-03/3D	
Nobel Biocare® Replace® Select	Ø3.5	M. Rot.Schutz	0.5	IPD/AC-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/AC-IN-02/3D	
		O. Rot.Schutz	0.5	IPD/AC-IN-01/3D	IPD/KT-IN-01
	Ø4.3	M. Rot.Schutz	1.5	IPD/AC-IN-03/3D	IPD/KT-IN-01
			0.5	IPD/AC-IR-00/3D	
			1.5	IPD/AC-IR-02/3D	
		O. Rot.Schutz	0.5	IPD/AC-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/AC-IR-03/3D	
			0.5	IPD/AC-IR-05/3D	
	Ø5.0	M. Rot.Schutz	0.5	IPD/AC-IW-00/3D	IPD/KT-IW-00
			1.5	IPD/AC-IW-02/3D	
		O. Rot.Schutz	0.5	IPD/AC-IW-01/3D	IPD/KT-IW-01
Nobel Biocare® Nobel Active®	Ø3.0	M. Rot.Schutz	0.8	IPD/AD-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/AD-IT-02/3D	
			2.5	IPD/AD-IT-04/3D	
		O. Rot.Schutz	0.8	IPD/AD-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/AD-IT-03/3D	
			2.5	IPD/AD-IT-05/3D	
	Ø3.5	M. Rot.Schutz	0.5	IPD/AD-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/AD-IN-02/3D	
			2.5	IPD/AD-IN-04/3D	
		O. Rot.Schutz	3.5	IPD/AD-IN-06/3D	IPD/KT-IN-01
			0.5	IPD/AD-IN-01/3D	
			1.5	IPD/AD-IN-03/3D	
Nobel Biocare® Nobel Active®	Ø4.3	M. Rot.Schutz	2.5	IPD/AD-IN-05/3D	IPD/KT-IN-01
			3.5	IPD/AD-IN-07/3D	
			0.5	IPD/AD-IR-00/3D	IPD/KT-IR-00
		O. Rot.Schutz	1.5	IPD/AD-IR-02/3D	
			2.5	IPD/AD-IR-04/3D	
			0.5	IPD/AD-IR-01/3D	IPD/KT-IR-01
	Ø4.3	O. Rot.Schutz	1.5	IPD/AD-IR-03/3D	
			2.5	IPD/AD-IR-05/3D	
			0.5	IPD/AD-IR-06/3D	
Biomet® 3i® External™ 3i®	Ø3.4	M. Rot.Schutz	1.5	IPD/BA-IN-02/3D	IPD/KT-IN-00
		O. Rot.Schutz		IPD/BA-IN-03/3D	IPD/KT-IN-01
	Ø5.0	M. Rot.Schutz	0.5	IPD/BA-IW-00/3D	IPD/KT-IW-00
			1.5	IPD/BA-IW-02/3D	
	Ø5.1	O. Rot.Schutz	0.5	IPD/BA-IW-01/3D	IPD/KT-IW-01
			1.5	IPD/BA-IW-03/3D	
Microdent® Microdent® Universal™	Ø5.1	O. Rot.Schutz	1.5	IPD/BA-IW-03/3D	IPD/KT-IW-01

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz			
Biomet® 3i® Certain®	Ø3.4	M. Rot.Schutz	0.5	IPD/BB-IN-00/3D	IPD/KT-IN-00			
			1.5	IPD/BB-IN-02/3D				
			2.5	IPD/BB-IN-04/3D				
		O. Rot.Schutz	0.5	IPD/BB-IN-01/3D	IPD/KT-IN-01			
			1.5	IPD/BB-IN-03/3D				
			2.5	IPD/BB-IN-05/3D				
	Ø4.1	M. Rot.Schutz	0.5	IPD/BB-IR-00/3D	IPD/KT-IR-00			
			1.5	IPD/BB-IR-02/3D				
			0.5	IPD/BB-IR-01/3D		IPD/KT-IR-01		
		1.5	IPD/BB-IR-03/3D					
		Ø5.0	M. Rot.Schutz	0.5	IPD/BB-IW-00/3D		IPD/KT-IW-00	
				1.5	IPD/BB-IW-02/3D			
	O. Rot.Schutz		0.5	IPD/BB-IW-01/3D	IPD/KT-IW-01			
	1.5	IPD/BB-IW-03/3D						
	Klockner® SK2-NK2	Ø4.2	M. Rot.Schutz	0.5	IPD/CA-IR-00/3D	IPD/KT-IW-00		
2.25				IPD/CA-IR-04/3D	IPD/KT-IR-00			
O. Rot.Schutz			0.5	IPD/CA-IR-01/3D	IPD/KT-IW-01			
			2.25	IPD/CA-IR-05/3D	IPD/KT-IR-01			
Ø4.5		M. Rot.Schutz	0.5	IPD/CB-IR-00/3D	IPD/KT-IR-00			
		O. Rot.Schutz		IPD/CB-IR-01/3D	IPD/KT-IR-01			
Klockner® Vega®	Ø3.5	M. Rot.Schutz	0.7	IPD/CD-IN-00/3D	IPD/KT-IN-00			
			1.5	IPD/CD-IN-02/3D				
			2.5	IPD/CD-IN-04/3D				
		O. Rot.Schutz	0.7	IPD/CD-IN-01/3D	IPD/KT-IN-01			
			1.5	IPD/CD-IN-03/3D				
			2.5	IPD/CD-IN-05/3D				
	Ø4.1/Ø4.5	M. Rot.Schutz	0.7	IPD/CD-IR-00/3D	IPD/KT-IR-00			
			1.5	IPD/CD-IR-02/3D				
			2.5	IPD/CD-IR-04/3D				
		O. Rot.Schutz	0.7	IPD/CD-IR-01/3D	IPD/KT-IR-01			
			1.5	IPD/CD-IR-03/3D				
			2.5	IPD/CD-IR-05/3D				
Straumann® Tissue Level	Ø4.8	M. Rot.Schutz	0.4	IPD/DA-IR-00/3D	IPD/KT-IR-00			
		Rotativa		IPD/DA-IR-01/3D	IPD/KT-IR-01			
	Ø6.5	M. Rot.Schutz		IPD/DA-IW-02/3D	IPD/KT-IW-00			
		O. Rot.Schutz		IPD/DA-IW-03/3D	IPD/KT-IW-01			
		Straumann® Bone Level		Ø3.3 NC	M. Rot.Schutz	0.5	IPD/DB-IN-00/3D	IPD/KT-IN-00
						1.5	IPD/DB-IN-02/3D	
2.5	IPD/DB-IN-04/3D							
O. Rot.Schutz	3.5		IPD/DB-IN-06/3D		IPD/KT-IN-01			
	0.5		IPD/DB-IN-01/3D					
	1.5		IPD/DB-IN-03/3D					
Ø4.1 RC	M. Rot.Schutz		2.5	IPD/DB-IN-05/3D	IPD/KT-IR-00			
			3.5	IPD/DB-IN-07/3D				
			0.5	IPD/DB-IR-00/3D		IPD/KT-IR-00		
	1.5		IPD/DB-IR-02/3D					
	2.5		IPD/DB-IR-04/3D					
	O. Rot.Schutz		3.5	IPD/DB-IR-06/3D	IPD/KT-IR-01			
0.5		IPD/DB-IR-01/3D						
1.5		IPD/DB-IR-03/3D						
2.5		IPD/DB-IR-05/3D						
			3.5	IPD/DB-IR-07/3D				

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Straumann® BLX®	RB/WB	M. Rot.Schutz	1.0	IPD/DC-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/DC-IR-02/3D	
			3.0	IPD/DC-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/DC-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/DC-IR-03/3D	
			3.0	IPD/DC-IR-05/3D	
Straumann® SRA®	Ø3.5	O. Rot.Schutz	0.2	IPD/DD-IN-01/3D	-
	Ø4.6			IPD/DD-IR-01/3D	
Astra® Osseospeed®	Ø3.0	M. Rot.Schutz	0.8	IPD/EA-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IN-02/3D	
			2.5	IPD/EA-IN-04/3D	
		O. Rot.Schutz	0.8	IPD/EA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/EA-IN-03/3D	
			2.5	IPD/EA-IN-05/3D	
	Ø3.5-4.0	M. Rot.Schutz	0.5	IPD/EA-IR-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IR-02/3D	
			2.5	IPD/EA-IR-04/3D	
		O. Rot.Schutz	3.5	IPD/EA-IR-06/3D	IPD/KT-IN-01
			0.5	IPD/EA-IR-01/3D	
			1.5	IPD/EA-IR-03/3D	
	Ø4.5-5.0	M. Rot.Schutz	2.5	IPD/EA-IR-05/3D	IPD/KT-IN-01
			3.5	IPD/EA-IR-07/3D	
			0.5	IPD/EA-IW-00/3D	IPD/KT-IR-00
		O. Rot.Schutz	1.5	IPD/EA-IW-02/3D	
			0.5	IPD/EA-IW-01/3D	
			1.5	IPD/EA-IW-03/3D	IPD/KT-IR-01
Astra® Evolution®	Ø3.6	M. Rot.Schutz	0.5	IPD/EB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EB-IN-02/3D	
			2.5	IPD/EB-IN-04/3D	
		O. Rot.Schutz	0.5	IPD/EB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/EB-IN-03/3D	
			2.5	IPD/EB-IN-05/3D	
Astra® Evolution®	Ø4.2	M. Rot.Schutz	0.5	IPD/EB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/EB-IR-02/3D	
			2.5	IPD/EB-IR-04/3D	
		O. Rot.Schutz	0.5	IPD/EB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/EB-IR-03/3D	
			2.5	IPD/EB-IR-05/3D	
Zimmer® Screw Vent®	Ø3.5	M. Rot.Schutz	0.5	IPD/FA-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/FA-IN-02/3D	
			2.5	IPD/FA-IN-04/3D	
		O. Rot.Schutz	0.5	IPD/FA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/FA-IN-03/3D	
			2.5	IPD/FA-IN-05/3D	
	Ø4.5	M. Rot.Schutz	0.5	IPD/FA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/FA-IR-02/3D	
			2.5	IPD/FA-IR-04/3D	
		O. Rot.Schutz	0.5	IPD/FA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/FA-IR-03/3D	
			2.5	IPD/FA-IR-05/3D	
	Ø5.7	M. Rot.Schutz	0.5	IPD/FA-IW-00/3D	IPD/KT-IW-00
		O. Rot.Schutz	0.5	IPD/FA-IW-01/3D	IPD/KT-IW-01

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Zimmer® Eztetic®	Ø3.1	M. Rot.Schutz	1.0	IPD/FC-IT-00/3D	IPD/KT-IN-00
			2.0	IPD/FC-IT-02/3D	
			3.0	IPD/FC-IT-04/3D	
		O. Rot.Schutz	1.0	IPD/FC-IT-01/3D	IPD/KT-IN-01
			2.0	IPD/FC-IT-03/3D	
			3.0	IPD/FC-IT-05/3D	
BTI® Tiny BTI®	Ø3.5	M. Rot.Schutz	1.5	IPD/GA-IN-02/3D	IPD/KT-IN-00
		O. Rot.Schutz		IPD/GA-IN-03/3D	IPD/KT-IN-01
BTI® Interna Universal®	Ø4.1	M. Rot.Schutz	0.5	IPD/GB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/GB-IR-02/3D	
			2.5	IPD/GB-IR-04/3D	
		O. Rot.Schutz	0.5	IPD/GB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/GB-IR-03/3D	
			2.5	IPD/GB-IR-05/3D	
BTI® Multi-im®	Ø4.1	O. Rot.Schutz	0.3	IPD/GC-IR-01/3D	-
Dentsply® Xive® Friadent®	Ø3.4	M. Rot.Schutz	1.5	IPD/IA-IN-02/3D	IPD/KT-IN-00
		O. Rot.Schutz	1.5	IPD/IA-IN-03/3D	IPD/KT-IN-01
	Ø3.8	M. Rot.Schutz	0.5	IPD/IA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/IA-IR-02/3D	
		O. Rot.Schutz	0.5	IPD/IA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/IA-IR-03/3D	
Dentsply® Ankylos®	C/X	M. Rot.Schutz	1.5	IPD/IB-IR-02/3D	IPD/KT-IR-00
			2.5	IPD/IB-IR-04/3D	
		O. Rot.Schutz	1.5	IPD/IB-IR-03/3D	IPD/KT-IR-01
			2.5	IPD/IB-IR-05/3D	
		M. Rot.Schutz	1.5	IPD/JA-IN-02/3D	IPD/KT-IN-00
			1.5	IPD/JA-IN-03/3D	
Camlog® Camlog®	Ø3.3	M. Rot.Schutz	1.5	IPD/JA-IR-00/3D	IPD/KT-IR-00
		O. Rot.Schutz	1.5	IPD/JA-IR-02/3D	
	Ø3.8	M. Rot.Schutz	0.5	IPD/JA-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/JA-IR-03/3D	
		O. Rot.Schutz	0.5	IPD/JA-IW-00/3D	IPD/KT-IR-00
			1.5	IPD/JA-IW-02/3D	
	Ø4.3	M. Rot.Schutz	0.5	IPD/JA-IW-01/3D	IPD/KT-IR-01
			1.5	IPD/JA-IW-03/3D	
		O. Rot.Schutz	0.5	IPD/JA-IX-00/3D	IPD/KT-IW-00
			1.5	IPD/JA-IX-02/3D	
	Ø5.0	M. Rot.Schutz	0.5	IPD/JA-IX-01/3D	IPD/KT-IW-01
		O. Rot.Schutz	1.5	IPD/JA-IX-03/3D	

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Camlog® Conelog®	Ø3.3	M. Rot.Schutz	0.5	IPD/JB-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/JB-IT-02/3D	
		O. Rot.Schutz	0.5	IPD/JB-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/JB-IT-03/3D	
	Ø3.8	M. Rot.Schutz	0.5	IPD/JB-IN-00/3D	IPD/KT-IR-00
			1.5	IPD/JB-IN-02/3D	
		O. Rot.Schutz	0.5	IPD/JB-IN-01/3D	IPD/KT-IR-01
			1.5	IPD/JB-IN-03/3D	
	Ø4.3	M. Rot.Schutz	0.5	IPD/JB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/JB-IR-02/3D	
		O. Rot.Schutz	0.5	IPD/JB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/JB-IR-03/3D	
BioHorizons® Tapered Internal	Ø3.0	M. Rot.Schutz	0.5	IPD/LB-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/LB-IT-02/3D	
		O. Rot.Schutz	2.5	IPD/LB-IT-04/3D	IPD/KT-IN-01
			0.5	IPD/LB-IT-01/3D	
	Ø3.5	M. Rot.Schutz	1.5	IPD/LB-IT-03/3D	IPD/KT-IN-00
			2.5	IPD/LB-IT-05/3D	
		O. Rot.Schutz	0.5	IPD/LB-IN-00/3D	IPD/KT-IN-01
			1.5	IPD/LB-IN-02/3D	
	Ø4.5	M. Rot.Schutz	0.5	IPD/LB-IN-04/3D	IPD/KT-IR-00
			1.5	IPD/LB-IN-01/3D	
		O. Rot.Schutz	2.5	IPD/LB-IN-03/3D	IPD/KT-IR-01
			0.5	IPD/LB-IR-00/3D	
	Ø5.7	M. Rot.Schutz	1.5	IPD/LB-IR-02/3D	IPD/KT-IR-00
			2.5	IPD/LB-IR-04/3D	
		O. Rot.Schutz	0.5	IPD/LB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/LB-IR-03/3D	
Sweden y Martina® Outlink®	Ø3.5	M. Rot.Schutz	2.5	IPD/LB-IR-05/3D	IPD/KT-IN-00
			0.5	IPD/LB-IN-05/3D	
		O. Rot.Schutz	1.5	IPD/LB-IW-00/3D	IPD/KT-IW-00
			0.5	IPD/LB-IW-01/3D	
	Ø5.0	M. Rot.Schutz	1.5	IPD/MA-IN-02/3D	IPD/KT-IN-00
			2.5	IPD/MA-IN-04/3D	
		O. Rot.Schutz	1.5	IPD/MA-IN-03/3D	IPD/KT-IN-01
			2.5	IPD/MA-IN-05/3D	
Sweden y Martina® Premium™ Kohnno®	Ø3.3	M. Rot.Schutz	0.5	IPD/BA-IW-00/3D	IPD/KT-IW-00
			1.5	IPD/BA-IW-02/3D	
		O. Rot.Schutz	0.5	IPD/BA-IW-01/3D	IPD/KT-IW-01
			1.5	IPD/BA-IW-03/3D	
	Ø3.8	M. Rot.Schutz	0.5	IPD/MB-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/MB-IT-02/3D	
		O. Rot.Schutz	0.5	IPD/MB-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/MB-IT-03/3D	
	Ø3.8	M. Rot.Schutz	0.7	IPD/MB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/MB-IN-02/3D	
		O. Rot.Schutz	0.7	IPD/MB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/MB-IN-03/3D	

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Sweden y Martina® Premium™ Kohno®	Ø4.25	M. Rot.Schutz	0.7	IPD/MB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/MB-IR-02/3D	
		O. Rot.Schutz	0.7	IPD/MB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/MB-IR-03/3D	
Phibo® TSH®	S2	M. Rot.Schutz	1.5	IPD/NA-IN-02/3D	IPD/KT-IN-00
		O. Rot.Schutz		IPD/NA-IN-03/3D	IPD/KT-IN-01
	S5	M. Rot.Schutz	1.3	IPD/NA-IW-02/3D	IPD/KT-IW-00
		O. Rot.Schutz		IPD/NA-IW-03/3D	IPD/KT-IW-01
Osstem® TSIII®	Ø3.5 Mini	M. Rot.Schutz	1.0	IPD/OB-IN-00/3D	IPD/KT-IN-00
			2.0	IPD/OB-IN-02/3D	
			3.0	IPD/OB-IN-04/3D	
		O. Rot.Schutz	1.0	IPD/OB-IN-01/3D	IPD/KT-IN-01
			2.0	IPD/OB-IN-03/3D	
	Ø4.0 Regular	M. Rot.Schutz	3.0	IPD/OB-IN-05/3D	IPD/KT-IN-01
			1.0	IPD/OB-IR-00/3D	
			2.0	IPD/OB-IR-02/3D	
		O. Rot.Schutz	3.0	IPD/OB-IR-04/3D	IPD/KT-IR-00
			1.0	IPD/OB-IR-01/3D	
Biotech® Dental Kontakt®	RP	M. Rot.Schutz	2.0	IPD/OB-IR-03/3D	IPD/KT-IR-01
			3.0	IPD/OB-IR-05/3D	
			1.0	IPD/QB-IR-00/3D	
		O. Rot.Schutz	2.0	IPD/QB-IR-02/3D	IPD/KT-IR-00
			3.0	IPD/QB-IR-04/3D	
			1.0	IPD/QB-IR-01/3D	
Neodent® Helix® HE	Ø3.3	M. Rot.Schutz	2.0	IPD/QB-IR-03/3D	IPD/KT-IR-01
			3.0	IPD/QB-IR-05/3D	
		O. Rot.Schutz	1.0	IPD/QB-IR-01/3D	-
	Ø5.0	M. Rot.Schutz	1.3	IPD/RA-IN-02/3D	-
			2.5	IPD/RA-IN-04/3D	
		O. Rot.Schutz	1.3	IPD/RA-IN-03/3D	-
Neodent® Grand Morse™ GM	GM	M. Rot.Schutz	2.5	IPD/RA-IN-05/3D	-
			1.3	IPD/RA-IW-02/3D	
			1.3	IPD/RA-IW-03/3D	
			1.3	IPD/RA-IW-03/3D	
		O. Rot.Schutz	1.0	IPD/RB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/RB-IR-02/3D	
			3.0	IPD/RB-IR-04/3D	
			4.0	IPD/RB-IR-06/3D	
Neodent® GM Abutment	RP	M. Rot.Schutz	1.0	IPD/RB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/RB-IR-03/3D	
		O. Rot.Schutz	3.0	IPD/RB-IR-05/3D	IPD/KT-IR-01
			4.0	IPD/RB-IR-07/3D	
Neodent® GM Micro Abutment	NP	M. Rot.Schutz	0.5	IPD/RC-IR-00/3D	IPD/KT-IW-00
			0.5	IPD/RD-IN-00/3D	
		O. Rot.Schutz	0.3	IPD/RD-IN-01/3D	IPD/KT-IN-01
			0.5	IPD/SB-IN-00/3D	
Bego® Semados®	Ø3.25/3.75	M. Rot.Schutz	1.5	IPD/SB-IN-02/3D	IPD/KT-IN-00
			2.5	IPD/SB-IN-04/3D	
			0.5	IPD/SB-IN-01/3D	
		O. Rot.Schutz	1.5	IPD/SB-IN-03/3D	IPD/KT-IN-01
			2.5	IPD/SB-IN-05/3D	
			2.5	IPD/SB-IN-05/3D	

Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Bego® Semados®	Ø4.1	M. Rot.Schutz	0.5	IPD/SB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/SB-IR-02/3D	
			2.5	IPD/SB-IR-04/3D	
		O. Rot.Schutz	0.5	IPD/SB-IR-01/3D	IPD/KT-IR-01
			1.5	IPD/SB-IR-03/3D	
			2.5	IPD/SB-IR-05/3D	
	Ø4.5	M. Rot.Schutz	0.5	IPD/SB-IW-00/3D	IPD/KT-IR-00
			1.5	IPD/SB-IW-02/3D	
			2.5	IPD/SB-IW-04/3D	
		O. Rot.Schutz	0.5	IPD/SB-IW-01/3D	IPD/KT-IR-01
			1.5	IPD/SB-IW-03/3D	
			2.5	IPD/SB-IW-05/3D	
MIS® Seven®	Narrow Ø3.3	M. Rot.Schutz	0.5	IPD/TA-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/TA-IT-02/3D	
			2.5	IPD/TA-IT-04/3D	
		O. Rot.Schutz	0.5	IPD/TA-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/TA-IT-03/3D	
			2.5	IPD/TA-IT-05/3D	
	Standard Ø3.75/4.2	M. Rot.Schutz	0.5	IPD/FA-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/FA-IN-02/3D	
			2.5	IPD/FA-IN-04/3D	
		O. Rot.Schutz	0.5	IPD/FA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/FA-IN-03/3D	
			2.5	IPD/FA-IN-05/3D	
Ancho Ø5.0/6.0	M. Rot.Schutz	0.5	IPD/FA-IR-00/3D	IPD/KT-IR-00	
		1.5	IPD/FA-IR-02/3D		
		2.5	IPD/FA-IR-04/3D		
	O. Rot.Schutz	0.5	IPD/FA-IR-01/3D	IPD/KT-IR-01	
		1.5	IPD/FA-IR-03/3D		
		2.5	IPD/FA-IR-05/3D		
MIS® C1/V3®	Standard	M. Rot.Schutz	0.7	IPD/TB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/TB-IR-02/3D	
			2.5	IPD/TB-IR-04/3D	
		O. Rot.Schutz	3.5	IPD/TB-IR-06/3D	IPD/KT-IR-01
			0.7	IPD/TB-IR-01/3D	
			1.5	IPD/TB-IR-03/3D	
DIO® UFII	NP	M. Rot.Schutz	2.5	IPD/TB-IR-05/3D	IPD/KT-IR-01
			3.5	IPD/TB-IR-07/3D	
			1.0	IPD/UB-IN-00/3D	
		2.0	IPD/UB-IN-02/3D		
		3.0	IPD/UB-IN-04/3D		
		O. Rot.Schutz	1.0	IPD/UB-IN-01/3D	IPD/KT-IN-01
2.0	IPD/UB-IN-03/3D				
3.0	IPD/UB-IN-05/3D				
DIO® UFII	RP	M. Rot.Schutz	1.0	IPD/UB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/UB-IR-02/3D	
			3.0	IPD/UB-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/UB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/UB-IR-03/3D	
			3.0	IPD/UB-IR-05/3D	



Kompatibilitäten	Plattform	Verbindung	GH (mm)	Artikelnr	Schneidaufsatz
Megagen® AnyOne®	RP	M. Rot.Schutz	1.0	IPD/WA-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/WA-IR-02/3D	
			3.0	IPD/WA-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/WA-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/WA-IR-03/3D	
			3.0	IPD/WA-IR-05/3D	
Megagen® AnyRidge®	RP	M. Rot.Schutz	1.0	IPD/WB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/WB-IR-02/3D	
			3.0	IPD/WB-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/WB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/WB-IR-03/3D	
			3.0	IPD/WB-IR-05/3D	
Dentium® Implantium® / Superline™	RP	M. Rot.Schutz	1.0	IPD/XA-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/XA-IR-02/3D	
			3.0	IPD/XA-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/XA-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/XA-IR-03/3D	
			3.0	IPD/XA-IR-05/3D	
Medentis® ICX	RP	M. Rot.Schutz	1.0	IPD/YB-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/YB-IR-02/3D	
			3.0	IPD/YB-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/YB-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/YB-IR-03/3D	
			3.0	IPD/YB-IR-05/3D	
Global D® In-Kone®	RP	M. Rot.Schutz	1.0	IPD/1A-IR-00/3D	IPD/KT-IR-00
			2.0	IPD/1A-IR-02/3D	
			3.0	IPD/1A-IR-04/3D	
		O. Rot.Schutz	1.0	IPD/1A-IR-01/3D	IPD/KT-IR-01
			2.0	IPD/1A-IR-03/3D	
			3.0	IPD/1A-IR-05/3D	
Anthogyr® Axiom® BL	RP	M. Rot.Schutz	1.5	IPD/2A-IR-02/3D	IPD/KT-IR-00
			2.5	IPD/2A-IR-04/3D	
			3.5	IPD/2A-IR-06/3D	
		O. Rot.Schutz	1.5	IPD/2A-IR-03/3D	IPD/KT-IR-01
			2.5	IPD/2A-IR-05/3D	
			3.5	IPD/2A-IR-07/3D	
Galimplant® Multi-Posición Estético®	Ø3.8	M. Rot.Schutz	0.2	IPD/3A-IR-00/3D	-
				IPD/3A-IR-02/3D	
		O. Rot.Schutz		IPD/3A-IR-01/3D	
				IPD/3A-IR-03/3D	

BE SIMPLY DIGITAL

IPD Scanbody



NT Scanbody



OUTPUT

preform Blank



Multi Unit



Digital Analog



Kronen & Brücken



Direkt Implantat



DIM Analog



CAD Lösungen für alle führenden Software-Systeme, kompatibel mit dem NT 3D Guide und dem IPD Scan Abutment System.

In Zusammenarbeit mit unserem Partner, der IPD Dental Group, bietet DENTALimpulse Ihnen präzise, sichere und preiswerte Lösungen für Ihre täglichen, digitalen Herausforderungen an.

Das Ziel von IPD ist es, Ihren digitalen Workflow zu vereinfachen und zu verbessern. Ein komplettes System, entwickelt und getestet, um Ihnen einen präzisen und sicheren digitalen Ablauf zu garantieren. Ob der 3D

Guide Scanbody von NT oder das innovative extra/intra-orale IPD Scan Abutment System:

Mit der ipd CAD-Bibliothek arbeiten Sie, im digitalen Workflow wie gewohnt, in Ihrer CAM Software (z.B. exocad, 3Shape, Blender etc.). So können Sie MDR-zertifizierte Titanbasen, Direktversorgungen auf Implantaten, Kronen und Brücken, sowie Scan Transfer preiswert und sicher anbieten.

